3870M Large Capacity Sterilizer







3870M

Item No	Chamber Dimensions Ø x L in(mm)	Volume	Overall Dimensions DxWxH in(mm)	Standard Cassettes	Tray Dimensions DxWxH in(mm)	No. of Trays	Standard Unwrapped Cycle Time	Voltage(V) Freq.(Hz)	Power(W) Current(A)	Shipping Weight Ibs(kg)
3870M with star	15"x30" d (384x760)	22gal/85L	34.5″x26″x20.6″ (875x665x525)	15 Full	26.6"x11.3"x1" (675x286x25) 26.6"x13.8"x1" (675x350x25)	2	Cold: 32 mins Hot: 23 mins	230V 50/60	3000W 13.6A	209 (95)

Features

- Automatic shut off at the end of both the sterilization and dry cycles
- · A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Drain valve is located on the front, allowing for quick and easy draining of water reservoir
- · Dual safety thermostat to protect against overheating
- · Supplied with a 3 shelf stand
- 1 year parts and labor warranty

Control Panel





MERCOFRAMES OPTICAL CORP

5555 Nw 74 Ave. Miami, Fl 33166 305-882-0120 www.mercoframes.com sales@mercoframes.net



3870EA Large Capacity Sterilizer







3870EA

Item No.	Chamber Dimensions Ø x L in(mm)	Volume	Overall Dimensions DxWxH in(mm)	Standard Cassettes	Tray Dimensions DxWxH in(mm)	No. of Trays	Standard Unwrapped Cycle Time	Voltage(V) Freq.(Hz)	Power(W) Current(A)	Shipping Weight Ibs(kg)
3870EA with stand	15″x30″ (384x760)	22gal/85L	34.5″x26″x20.6″ (875x665x525)	15 Full	26.6"x11.3"x1" (675x286x25) 26.6"x13.8"x1" (675x350x25)	2	Cold: 32 mins Hot: 23 mins	230V 50/60	3000W 13.6A	213 (97)

Features

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- A long life electro polished chamber and door
- New brighter and easier to read panel
- Program parameters can be customized and stored
- · Supplied with a 3 shelf stand
- 2 year parts and labor warranty

Control Panel



- 1. Program Selection
- 2. Detailed Digital Display
- Programmable Keys
- 4. System Indicators

Optional Internal Printer



Documents: date, time, temperature, pressure.

Control of the contro

Clean & Simple™ Ultrasonic Cleaners and Tablets



Instrument Pouch Rack



Steam Distiller



Chamber Brite Autoclave Cleaner